



Recombinant Human Protein WFDC9 (WFDC9)

Product Code	CSB-MP822776HU
Abbreviation	WFDC9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEX5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SFWNKDP FLDMIRETEQ CWVQPPYKYC EKRCTKIMTC VRPNHTCCWT YCGNICLDNE EPLKSMLNP
Source	Mammalian cell
Target Names	WFDC9
Protein Names	Recommended name: Protein WFDC9
Expression Region	24-89
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The WAP-type four-disulfide core (WFDC) domain, or WAP signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many members of the WFDC domain family. This gene encodes a protein which contains a WFDC domain, and is thus a member of the WFDC domain family. This gene and several other gene family members are clustered at 20q13.12.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.