



# Recombinant Human TNFAIP3-interacting protein 2 (TNIP2)

<b>Product Code</b>	CSB-EP822786HU-B
<b>Abbreviation</b>	TNIP2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NFZ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSRDPGSGGW EEAPRAAAAL CTLYHEAGQR LRRLQDQLAA RDALIARLRA RLAALEGDAA PSLVDALLEQ VARFREQLRQ QEGGAAEAQM RQEIERLTER LEEKEREMQQ LLSQPQHERE KEVLLRRSM AEGERARAAS DVLCRSLANE THQLRRTLTA TAHMCQHAK CLDERQHAQR NVGERSPDQS EHTDGHTSVQ SVIEKLQEEEN RLLKQKVTHV EDLNAKWQRY NASRDEYVRG LHAQLRGLQI PHEPELMRKE ISRLNRQLEE KINDCAEVKQ ELAASRTARD AALERVQMLE QQILAYKDDF MSERADRERA QSRIQELEEK VASLLHQVSW RQDSREPDAG RIHAGSKTAK YLAADALELM VPGGWRPGTG SQQPEPPAEG GHPGAAQRGQ GDLQCPHCLQ CFSDEQGEEL LRHVAECCQ
<b>Source</b>	E.coli
<b>Target Names</b>	TNIP2
<b>Protein Names</b>	Recommended name: TNFAIP3-interacting protein 2 Alternative name(s): A20-binding inhibitor of NF-kappa-B activation 2 Short name= ABIN-2 Fetal liver LKB1-interacting protein
<b>Expression Region</b>	1-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.