



Recombinant Human TNFAIP3-interacting protein 2 (TNIP2)

Product Code	CSB-MP822786HU
Abbreviation	TNIP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NFZ5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSRDPGSGGW EEAPRAAAL CTLYHEAGQR LRRLQDQLAA RDALIARLRA RLAALEGDAA PSLVDALLEQ VARFREQLRQ QEGGAAEAQM RQEIERLTER LEEKEREMQQ LLSQPQHERE KEVLLRRSM AEGERARAAS DVLCRSLANE THQLRRTLTA TAHMCQHAK CLDERQHAQR NVGERSPDQS EHTDGHTSVQ SVIEKLQEEEN RLLKQKVTHV EDLNAKWQRY NASRDEYVRG LHAQLRGLQI PHEPELMRKE ISRLNRQLEE KINDCAEVKQ ELAASRTARD AALERVQMLE QQILAYKDDF MSERADRERA QSRIQELEEK VASLLHQVSW RQDSREPDAG RIHAGSKTAK YLAADALELM VPGGWRPGTG SQQPEPPAEG GHPGAAQRGQ GDLQCPHCLQ CFSDEQGEEL LRHVAECCQ
Source	Mammalian cell
Target Names	TNIP2
Protein Names	Recommended name: TNFAIP3-interacting protein 2 Alternative name(s): A20-binding inhibitor of NF-kappa-B activation 2 Short name= ABIN-2 Fetal liver LKB1-interacting protein
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.