



Recombinant Human Tubulin--tyrosine ligase (TTL)

Product Code	CSB-MP822788HU
Abbreviation	TTL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NG68
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MYTFVVRDEN SSVYAEVSRL LLATGHWKRL RRDNPRFNLM LGERNRLPFG RLGHEPGLVQ LVNYYRGADK LCRKASLVKL IKTPELAES CTWFPESYVI YPTNLKTPVA PAQNGIQPPI SNSRTEREF FLASYNRKKE DGEGNVWIAK SSAGAKGEGI LISSEASELL DFIDNQGQVH VIQKYLEHPL LLEPGHRKFD IRSWVLVDHQ YNIYLYREGV LRTASEPYHV DNFQDKTCHL TNHCIQKEYS KNYGKYEEN EMFFKEFNQY LTSALNITLE SSILLQIKHI IRNCLLSVEP AISTKHLPYQ SFQLFGDFM VDEELKVWLI EVNGAPACAQ KLYAELCQGI VDIAISSVFP PPDVEQPQTQ PAAFIKL
Source	Mammalian cell
Target Names	TTL
Protein Names	Recommended name: Tubulin--tyrosine ligase Short name= TTL EC= 6.3.2.25
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.