



Recombinant Human Protein NEDD1 (NEDD1)

Product Code	CSB-EP822820HU-B
Abbreviation	NEDD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NHV4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQENLRFASS GDDIKIWDAS SMTLVDFKFNPT HTSPHGISSI CWSSNNNFFLV TASSSGDKIV VSSCKCKPVP LLELAEGQKQ TCVNLNSTSM YLVSGGLNNT VNIWDLKSKR VHRSLKDHKD QVTCVTYNWN DCYIASGSLG GEIILHSVTT NLSSTPFGHG SNQSVRHLKY SLFKKSLGVS VSDNGIVTLW DVNSQSPYHN FDSVHKAPAS GICFSPVNEL LFVTIGLDRK IILYDTSSKK LVKTLVADTP LTAVDFMPDG ATLAIGSSRG KIYQYDLRML KSPVKTISAH KTSVQCFIAFQ YSTVLTKSSL NKGCSNKPTT VNKRSVNVNA ASGGVQNSGI VREAPATSIA TVLPQPM TSA MGKGTAVVQE KAGLPRSINT DTLSKETDSG KNQDFSSFDD TGKSSLGDMF SPIRDDAVVN KGSDSIGKDG DGFDFLPQLN SVFPPRKNPV TSSTSVLHSS PLNVFMGSPG KEENENRDLT AESKKIYMGK QESKDSFKQL AKLVTSGAES GNLNTSPSSN QTRNSEKFEK PENEIEAQLI CEPPINGSST PNPKIASSVT AGVASSLSEK IADSIGNNRQ NAPLTSIQIR FIQNMIQETL DDFREACHRD IVNLQVEMIK QFHMQLNEMH SLLERYSVNE GLVAEIERLR EENKRLRAHF</p>
Source	E.coli
Target Names	NEDD1
Protein Names	Recommended name: Protein NEDD1 Alternative name(s): Neural precursor cell expressed developmentally down-regulated protein 1 Short name= NEDD-1
Expression Region	1-660
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.