



Recombinant *Xenopus laevis* Mitotic spindle assembly checkpoint protein MAD2B (mad2l2)

Product Code	CSB-EP823126XBE
Abbreviation	mad2l2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8QFR4
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MTTLTRQDLN FGQVVADILC EFLEVAVHLI LYVREVYPTG IFQKRKKYNV PVQMSCHPEL NRYIQDTLHC VKPLIEKNDV EKVVVVILDK EHHPVERFVF EIAQPPLLSI SSDSLLSHVE QLLRAFILKI SVCDAVLNPPGCTFTLLV HTREAATRNM EKIQVIKDFP WILADEQDVH MQEPRLIPLK TMTSDILKMQ LYVEERAQKS T
Source	E.coli
Target Names	mad2l2
Protein Names	Recommended name: Mitotic spindle assembly checkpoint protein MAD2B Alternative name(s): Mad2l2-A protein Mitotic arrest defective protein 2B Mitotic arrest deficient 2-like protein 2 Short name= MAD2-like protein 2 Short name=
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.