



Recombinant Hordeum vulgare Chalcone--flavonone isomerase (CHI)

Product Code	CSB-MP823329HWQ
Abbreviation	CHI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8S3X0
Product Type	Recombinant Protein
Immunogen Species	Hordeum vulgare (Barley)
Purity	>85% (SDS-PAGE)
Sequence	MAVSELEVDG VVFPPLARPP GSAHAHFLAG AGVRGMEIGG HFIKFTAIGV YLQADAAVSA LAAKWAGKPA ADLASDAFF RDVVTGEFEK FTRVTMILPL TGAQYSDKVT ENCVAYWKAA GVYTDAAEAAA VDKFKEAFGP HSFAPGASIL FTHSPAGVLT VAFSKDSSVP ESGGVAIENA RLCEAVLESI IGEHGVSPAA KLSLANRVAE LLKGA AHAGG EPAAEPVPS V
Source	Mammalian cell
Target Names	CHI
Protein Names	Recommended name: Chalcone--flavonone isomerase Short name= Chalcone isomerase EC= 5.5.1.6
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.