



Recombinant *Oryza sativa* subsp. japonica Ubiquitin-fold modifier-conjugating enzyme 1 (Os10g0205200, LOC_Os10g13800)

Product Code	CSB-YP823331OFG
Abbreviation	Os10g0205200, LOC_Os10g13800
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8S625
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEGWDKGTKS VVGEIPLLST RAGPRDGEAW RQRLKEEYRA LIAYTSVNKS KDNDWFRISA ANPEGTRWEG TCWYVHNLRR YEFPLQFDIP VAYPQVAPEI ELPTLDGKTH KMYRGGKICL TVHFKPLWAK NCPRFIAHA LCLGLAPWLA AEVPILVDSG MVKHKDDEAA PADAAAAASG SAAAS
Source	Yeast
Target Names	Os10g0205200
Protein Names	Recommended name: Ubiquitin-fold modifier-conjugating enzyme 1 Alternative name(s): Ufm1-conjugating enzyme 1
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.