



Recombinant *Macaca fascicularis* CD59 glycoprotein (CD59)

| | |
|--------------------------|---|
| Product Code | CSB-MP823365MOV |
| Abbreviation | CD59 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8SQ46 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey) |
| Purity | >85% (SDS-PAGE) |
| Sequence | LQCYN CPNPTTDCCKT AINCSSGFDTC LIARAGLQV YNQCWKFANCC NYNDISTLLK ESELYFCCK KDLCNFNEQL ES |
| Source | Mammalian cell |
| Target Names | CD59 |
| Protein Names | Recommended name: CD59 glycoprotein Alternative name(s): MAC-inhibitory protein Short name= MAC-IP Membrane attack complex inhibition factor Short name= MACIF Protectin CD_antigen= CD59 |
| Expression Region | 26-102 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |