



# Recombinant Human Kelch domain-containing protein 4 (KLHDC4)

<b>Product Code</b>	CSB-BP823448HU
<b>Abbreviation</b>	KLHDC4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TBB5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGKKGKKEKK GRGAEKTAAK MEKKVSKRSR KEEEDLEALI AHFQTLDAKR TQTVELPCPP PSPRLNASLS VHPEKDELIL FGGEYFNGQK TFLYNELYVY NTRKDTWTKV DIPSPPPRRC AHQAVVVPQG GGQLWVFGGE FASPNGEQFY HYKDLWVLHL ATKTWEQVKS TGGPSGRSGH RMVAWKRQLI LFGGFHESTR DYIYYNDVYA FNLDTFTWSK LSPSGTGPTP RSGCQMSVTP QGGIVVYGGY SKQRVKKDVD KGTRHSDMFL LKPEDGREDK WWWTRMNSPG VKPTPRSGFS VAMAPNHQTL FFGGVCDEEE EESLSGEFFN DLYFYDATRN RWFEGQLKGP KSEKKRRRRG RKEEPEGGSR PACGGAGTQG PVQLVKEVVA EDGTVVTIKQ VLTAPGSAGQ PRSEDEDSLE EAGSPAPGPC PRSNAMLAVK HGVLYVYGGM FEAGDRQVTL SDLHCLDLHR MEAWKALVEM DPETQEWLEE TDSEEDSEEV EGAEGGVDDE DSGEESGAED
<b>Source</b>	Baculovirus
<b>Target Names</b>	KLHDC4
<b>Protein Names</b>	Recommended name: Kelch domain-containing protein 4
<b>Expression Region</b>	1-520
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.