



Recombinant Human WD repeat-containing protein 20 (WDR20)

Product Code	CSB-MP823451HU
Abbreviation	WDR20
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TBZ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ATEGGGKEM NEIKTQFTTR EGLYKLLPHS EYSRPNRVPF NSQGSNPVRV SFVNLNDQSG NGDRLCFNVG RELYFYIYKG VRKAADLSKP IDKRIYKGTQ PTCHDFNHLT ATAESVSLLV GFSAGQVQLI DPIKKETSKL FNEERLIDKS RVTCVKWVPG SESLFLVAHS SGNMYLYNVE HTC GTTAPHY QLLKQGESFA VHTCKSKSTR NPLLKWTVGE GALNEFAFSP DGKFLACVSQ DGFLRVFNFD SVELHGTMKY YFGGLLCVCW SPDGKYIVTG GEDDLVTWVS FVDCRVIARG HGHKSWVSVV AFDPYTTSVE EGDPMEFSGS DEDFQDLLHF GRDRANSTQS RLSKRNSTDS RPSVSTYRFG SVGQDTQLCL WDLTEDILFP HQPLSRARTH TNVMNATSP AGSNGNSVTT PGNSVPPPLP RSNLPHSAV SNAGSKSSVM DGAIASGVSK FATLSLHDKR ERHHEKDKHR NHSMGHISSK SSDKLNLVTK TKTDPAKTLG TPLCPMEDV PLLEPLICKK IAHERLTVLI FLEDCIVTAC QEGFICTWGR PGKVVSFNP
Source	Mammalian cell
Target Names	WDR20
Protein Names	Recommended name: WD repeat-containing protein 20 Alternative name(s): Protein DMR
Expression Region	2-569
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a WD repeat-containing protein. Multiple alternatively spliced transcript variants have been found for this gene, but the full-length nature of some variants has not been defined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.