



Recombinant *Lapemis hardwickii* 60S ribosomal protein L36 (RPL36)

Product Code	CSB-YP823710LLG
Abbreviation	RPL36
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8UW19
Product Type	Recombinant Protein
Immunogen Species	<i>Hydrophis hardwickii</i> (Hardwick's spine-bellied seasnake) (<i>Lapemis hardwickii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MAIRYPMAVG LNKGHKVTKN VSKPRHCRRR GRLTKHTKFV RDMIREVCGF APYERRAMEL LKVSCKDKRAL KFIKKRVGTH IRAKRKREEL SNVLAAMRKA AAKKD
Source	Yeast
Target Names	RPL36
Protein Names	Recommended name: 60S ribosomal protein L36
Expression Region	1-105
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.