



Recombinant Monkeypox virus Intermediate transcription factor 3 large subunit (VITF3L)

Product Code	CSB-EP823719MHV-B
Abbreviation	VITF3L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8V4V4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Monkeypox virus
Purity	>85% (SDS-PAGE)
Sequence	MDNLFTFLHE IEDRYARTIF NFHLISCDEI GDIYGLMKER ISSEDMFDNI VYNKDIHPAI KKLVYCDIQL TKHIINQNTY PVFNDSSQVK CCHYFDINSD NSNISSRTVE IFESEKSSLV SYIKTTNKKR KVNYGEIKKT VHGGTNANYF SGKKSDEYLS TTVRSNINQP WIKTISKRMR VDIINHSIVT RGKSSILQTI EIIFTNRTCV KIFKDSTMHI ILSKDKDEKG CINMIDKLFY VYYNLFLLFE DIIQNEYFKE VANVVNHVLM ATALDEKLFL IKKMAEHDVY GVSNFKIGMF NLTFIKSLDH TVFPSLLDED SKIKFFKGKK LNIVALRSLE DCTNYVTKSE NMIEMMKERS TILNSIDIET ESVDRLKELL LK
Source	E.coli
Target Names	VITF3L
Protein Names	Recommended name: Intermediate transcription factor 3 large subunit Alternative name(s): VITF-3 45 kDa subunit
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.