



# Recombinant Arabidopsis thaliana UDP-glycosyltransferase 73B1 (UGT73B1)

<b>Product Code</b>	CSB-EP823813DOA
<b>Abbreviation</b>	UGT73B1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8VZE9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGTPVEVSKL HLLFPFMAH GHMIPTLDMA KLFATKGAKS TILTTPLNAK LFFEKPIKSF NQDNPGLEDI TIQILNFPCT ELGLPDGCEN TDFIFSTPDL NVGDLSQKFL LAMKYFEEPL EELLVTMRPD CLVGNMFFPW STKVAEKFGV PRLVFHGTGY FSLCASHCIR LPKNVATSSE PFVIPDLPGD ILITEEQVME TEESVMGRF MKAIRDSEED SFGVLVNSFY ELEQAYSDFY KSFVAKRAWH IGPLSLGNRK FEEKAERGKK ASIDEHECLK WLDSKKCDSV IYMAFGTMSS FKNEQLIEIA AGLDMSGHDF VVVVNRKGSQ VEKEDWLPEG FEEKTKGKGL IIRGWAPQVL ILEHKAIGGF LTHCGWNSLL EGVAAGLPMV TWPVGAEQFY NEKLVTQVLK TGVSVGVKMM MQVVGDFISR EKVEGAVREV MVGEERRKRA KELAEMAKNA VKEGGSSDLE VDRLMEELTL VKLQKEKV
<b>Source</b>	E.coli
<b>Target Names</b>	UGT73B1
<b>Protein Names</b>	Recommended name: UDP-glycosyltransferase 73B1 EC= 2.4.1.-
<b>Expression Region</b>	1-488
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.