



Recombinant Arabidopsis thaliana GATA transcription factor 10 (GATA10)

Product Code	CSB-MP823815DOA
Abbreviation	GATA10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VZP4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNWLPEAEAE EHLKGILSGD FFDGLTNHLD CPLEDIDSTN GEGDWVARFQ DLEPPPLDMF PALPSDLTSC PKGAARVRIP NNMIPALKQS CSSEALSGIN STPHQSSAPP DIKVSYLFQS LTPVSVLENS YGSLSTQNSG SQRLAFPVKG MRSKRRRPTT VRLSYLFPFE PRKSTPGESV TEGYYSSEQH AKKKRKIHLI THTESSTLES SKSDGIVRIC THCETITTPQ WRQGPSGPKT LCNACGVRFK SGRLVPEYRP ASSPTFIPSV HSNSHRKIIE MRKKDDEFDT SMIRSDIQKV KQGRKKMV
Source	Mammalian cell
Target Names	GATA10
Protein Names	Recommended name: GATA transcription factor 10
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.