



Recombinant Human Ras-related protein Rab-2B (RAB2B)

Product Code	CSB-BP823880HU
Abbreviation	RAB2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUD1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTYAYLFKYI IIGDTGVGKS CLLLQFTDKR FQPVHDLTIG VEFGARMVNI DGKQIKLQIW DTAGQESFRS ITRSYYRGAA GALLVYDITR RETFNHLTSW LEDARQHSSS NMVIMLIGNK SDLESRRDVK REEGEAFARE HGLIFMETS A KTACNVEEAF INTAKEIYRK IQQGLFDVHN EANGIKIGPQ QSISTSVGPS ASQRNSRDIG SNSGCC
Source	Baculovirus
Target Names	RAB2B
Protein Names	Recommended name: Ras-related protein Rab-2B
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.