



Recombinant Human AN1-type zinc finger protein 2B (ZFAND2B)

Product Code	CSB-MP823888HU
Abbreviation	ZFAND2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WV99
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEFPDLGAHC SEPSCQRLDF LPLKCDACSG IFCADHVAYA QHHCGSAYQK DIQVPVCPLC NVPVPVARGE PPDRAVGEHI DRDCRSDPAQ QKRKIFTNKC ERAGCRQREM MKLTCERCSSR NFCIKHRHPL DHDCSGEGHP TSRAGLAAIS RAQAVASTST VPSPSQTMPS CTSPSRATTR SPSWTAPPVI ALQNGLSEDE ALQRALEMSL AETKPQVPSC QEEEDLALAQ ALSASEAEYQ RQQAQSRSSK PSNCSLC
Source	Mammalian cell
Target Names	ZFAND2B
Protein Names	Recommended name: AN1-type zinc finger protein 2B Alternative name(s): Arsenite-inducible RNA-associated protein-like protein Short name= AIRAP-like protein
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.