



Recombinant Human Spindle and kinetochore-associated protein 2 (SKA2)

Product Code	CSB-MP823890HU
Abbreviation	SKA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WVK7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEAEVDKLEL MFQKAESDL D YIQYRLEYEI KTNHPDSASE KNPVTLKEL SVIKSRYQTL YARFKPVAVE QKESKSRICA TVKKTMMMIQ KLQKQTDLEL SPLTKEEKTA AEQFKFHMPD L
Source	Mammalian cell
Target Names	SKA2
Protein Names	Recommended name: Spindle and kinetochore-associated protein 2 Alternative name(s): Protein FAM33A
Expression Region	1-121
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.