



Recombinant Human Kelch repeat and BTB domain-containing protein 7 (KBTBD7)

Product Code	CSB-EP823893HU-B
Abbreviation	KBTBD7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WVZ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQSREDVPRS RRLASPRGGR RPKRISKPSV SAFFTGPEEL KDTAHS AALL AQLKSFYDAR LLC DVTIEVV TPGSGPGTGR LFSCNRNVLA AACPYFKSMF TGGMYESQQA SVTMHDVDAE SFEVLVDYCY TGRVSLSEAN VQRLYAASDM LQLEYVREAC ASFLARRLDL TNCTAILKFA DAFDHHKLR SQAQSYIAHNF KQLSRMGSIR EETLADTLA QLLAVLR LDS LDIESERTVC HVAVQWLEAA AKERGPSAAE VFKCVRWMHF TEEDQDYLEG LLTKPIVKKY CLDVIEGALQ MRYGDLLYKS LVPV PNSSSS SSSNSLVSA AENPPQRLGM CAKEMVIFFG HPRDPFLCYD PYSGDIYTMP SPLTSFAHTK TVTSSAVCVS PDHDIYLA AQ PRKDLWVYKP AQNSWQQLAD RLLCREGMDV AYLNGYIYIL GGRDPITGVK LKEVECYSVQ RNQWALVAPV PHSFYSELI VVQNYLYAVN SKRMLCYDPS HNMWLN CASL KRSDFQEACV FNDEIYCICD IPVMKVYNPA RGEWRRISNI PLDSETHNYQ IVNHDQKLLL ITSTTPQWKK NRVTVEYEDT REDQWINIGT MLG LLQFDSG FICLCARVYP SCLEPGQSFI TEEDDARSE S STEWDL DGFS ELDSESGSSS SFS DDEVVWQ VAPQRNAQDQ QGSL</p>
Source	E.coli
Target Names	KBTBD7
Protein Names	Recommended name: Kelch repeat and BTB domain-containing protein 7
Expression Region	1-684
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.