



# Recombinant Human Kelch repeat and BTB domain-containing protein 7 (KBTBD7)

<b>Product Code</b>	CSB-MP823893HU
<b>Abbreviation</b>	KBTBD7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WVZ9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>MQSREDVPRS RRLASPRGGR RPKRISKPSV SAFFTGPEEL KDTAHS AALL  AQLKSFYDAR LLC DVTIEVV TPGSGPGTGR LFSCNRNVLA AACPYFKSMF  TGGMYESQQA SVTMHDVDAE SFEVLVDYCY TGRVSLSEAN VQRLYAASDM  LQLEYVREAC ASFLARRLDL TNCTAILKFA DAFDHHKLR SQAQSYIAHNF  KQLSRMGSIR EETLADTLA QLLAVLR LDS LDIESERTVC HVAVQWLEAA  AKERGPSAAE VFKCVRWMHF TEEDQDYLEG LLTKPIVKKY CLDVIEGALQ  MRYGDLLYKS LVPV PNSSSS SSSNSLVSA AENPPQRLGM CAKEMVIFFG  HPRDPFLCYD PYSGDIYTMP SPLTSFAHTK TVTSSAVCVS PDHDIYLA AQ  PRKDLWVYKP AQNSWQQLAD RLLCREGMDV AYLNGYIYIL GGRDPITGVK  LKEVECYSVQ RNQWALVAPV PHSFYSELI VVQNYLYAVN SKRMLCYDPS  HNMWLN CASL KRSDFQEACV FNDEIYCICD IPVMKVYNPA RGEWRRISNI  PLDSETHNYQ IVNHDQKLLL ITSTTPQWKK NRVTVEYEDT REDQWINIGT  MLG LLQFDSG FICLCARVYP SCLEPGQSFI TEEDDARSES STEWDL DGFS  ELDSESGSSS SFS DDEVVWQ VAPQRNAQDQ QGSL</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KBTBD7
<b>Protein Names</b>	Recommended name: Kelch repeat and BTB domain-containing protein 7
<b>Expression Region</b>	1-684
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.