



Recombinant Human Transcriptional repressor p66-beta (GATAD2B)

Product Code	CSB-MP823905HU
Abbreviation	GATAD2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WXI9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDRMTEDALR LNLLKRSLDP ADERDDVLAK RLKMEGHEAM ERLKMLALLK RKDLANLEVP HELPTKQDGS GVKGYEEKLN GNLRPHGDNR TAGRPGKENI NDEPVDMSAR RSEPERGRLT PSPDIIVLS NEASSPRSSS RMEERLKAAN LEMFKGGKIE ERQQLIKQLR DELRLEEALR VLLKCLRQSQ LQKENVVQKT PVVQNAASIV QPSPAHVGGQ GLSKLPSRPG AQQVEPQNL TLQGHSVIRS ATNTTLPHML MSQRVIAPNP AQLQGQRGPP KPGLVRTTTP NMNPAINYQP QSSSSVPCQR TTSSAIYMNL ASHIQPGTVN RVSSPLPSPS AMTDAANSQA AAKLALRKQL EKTLLIPPP KPPAPLLHFL PSAANSEFIY MVGLEEVVQS VIDSQGKSCA SLLRVEPFVC AQCRTDFTPH WKQEKNKIL CEQCMTSNQK KALKAHTNR LKNAFVKALQ QEQEIEQLR QQAALSPPTA PAVSSVSKQE TIMRHHTLRQ APQPQSSLQR GIPTSARSML SNFAQAPQLS VPGGLLGMPG VNIAYLNTGI GGHKGPSLAD RQREYLLDMI PPRSISQSIG GQK
Source	Mammalian cell
Target Names	GATAD2B
Protein Names	Recommended name: Transcriptional repressor p66-beta Alternative name(s): GATA zinc finger domain-containing protein 2B p66/p68
Expression Region	1-593
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.