



# Recombinant Human Transcriptional repressor p66-beta (GATAD2B)

<b>Product Code</b>	CSB-YP823905HU
<b>Abbreviation</b>	GATAD2B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WXI9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDRMTEDALR LNLLKRSLDP ADERDDVLAK RLKMEGHEAM ERLKMLALLK RKDLANLEVP HELPTKQDGS GVKGYEEKLN GNLRPHGDNR TAGRPGKENI NDEPVDMSAR RSEPERGRLT PSPDIIVLS NEASSPRSSS RMEERLKAAN LEMFKGGKIE ERQQLIKQLR DELRLEEERL VLLKCLRQSQ LQKENVVQKT PVVQNAASIV QPSPAHVGQQ GLSKLPSRPG AQQVEPQNL TLQGHSVIRS ATNTTLPHML MSQRVIAPNP AQLQGQRGPP KPGLVRTTTP NMNPAINYQP QSSSSVPCQR TTSSAIYMNL ASHIQPGTVN RVSSPLPSPS AMTDAANSQA AAKLALRKQL EKTLLIPPP KPPAPLLHFL PSAANSEFIY MVGLEEVVQS VIDSQGKSCA SLLRVEPFVC AQCRTDFTPH WKQEKNKIL CEQCMTSNQK KALKAHTNR LKNAFVKALQ QEQEIEQLQ QQAALSPPTA PAVSSVSKQE TIMRHHTLRQ APQPQSSLQR GIPTSARSML SNFAQAPQLS VPGGLLGMPG VNIAYLNTGI GGHKGPSLAD RQREYLLDMI PPRSISQSIG GQK
<b>Source</b>	Yeast
<b>Target Names</b>	GATAD2B
<b>Protein Names</b>	Recommended name: Transcriptional repressor p66-beta Alternative name(s): GATA zinc finger domain-containing protein 2B p66/p68
<b>Expression Region</b>	1-593
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.