



Recombinant Human Intraflagellar transport protein 81 homolog (IFT81)

Product Code	CSB-EP823910HU
Abbreviation	IFT81
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WYA0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SDQIKFIMD SLNKEPFRKN YNLITFDSLE PMQLLQVLSL VLAEIDPKQL VDIREEMPEQ TAKRMLSLG ILKYKPSGNA TDMSTFRQGL VIGSKPVIYP VLHWLLQRTN ELKKRAYLAR FLIKLEVPSE FLQDETVA DT NKQYEELMEA FKTLHKEYEQ LKISGFSTAE IRKDISAMEE EKDQLIKRVE HLKKRVETAQ NHQWMLKIAR QLRVEKEREE YLAQQKQEQK NQLFHAVQRL QRVQNQLKSM RQAAADAKPE SLMKRLEEEI KFNLYMVTEK FPKELNKKK ELHFLQKVVS EPAMGHSDLL ELESKINEIN TEINQLIEKK MMRNEPIEGK LSLYRQQASI ISRKKEAKAE ELQEAKEKLA SLEREASVKR NQTREFDGTE VLKGDEFKRY VNKLRKSTV FKKKHQIIAE LKA EFGLLQR TEELLKQRHE NIQQQLQTME EKKGISGYSY TQEELERVSA LKSEVDEMKG RTLDDMSEMV KKLYSLVSEK KSALASVIKE LRQLRQKYQE LTQECDEKKS QYDSCAAGLE SNRSKLEQEV RRLREECLQE ESRYHYTNM IKNLEVQLRR ATDEM KAYIS SDQQEKRKAI REQYTKNTAE QENLGKKLRE KQKVIRESHG PNMKQAKMWR DLEQLMECKK QCFLKQQSQT SIGQVIQEGG EDRLIL
Source	E.coli
Target Names	IFT81
Protein Names	Recommended name: Intraflagellar transport protein 81 homolog Alternative name(s): Carnitine deficiency-associated protein expressed in ventricle 1 Short name= CDV-1
Expression Region	2-676
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.