



Recombinant Donkey Cytochrome c (CYCS)

| | |
|--------------------------|--|
| Product Code | CSB-BP006328DK |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P68097 |
| Product Type | Recombinant Protein |
| Immunogen Species | Equus asinus (Donkey) (Equus africanus asinus) |
| Purity | >85% (SDS-PAGE) |
| Sequence | GDVEKGKKI FVQKCAQCHT VEKGGKHKGTG PNLHGLFGRK TGQAPGFSYT DANKNKGITW KEETLMEYLE NPKKIYIPGTK MIFAGIKKKT EREDLIAYLK KATNE |
| Source | Baculovirus |
| Target Names | CYCS |
| Protein Names | Recommended name: Cytochrome c |
| Expression Region | 2-105 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | <p>This gene encodes cytochrome c, a component of the electron transport chain in mitochondria. The heme group of cytochrome c accepts electrons from the b-c1 complex and transfers electrons to the cytochrome oxidase complex.</p> <p>Cytochrome c is also involved in initiation of apoptosis. Upon release of cytochrome c to the cytoplasm, the protein binds apoptotic protease activating factor which activates the apoptotic initiator procaspase 9. Many cytochrome c pseudogenes exist, scattered throughout the human genome.</p> |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |