



# Recombinant Human Mitogen-activated protein kinase 3 (MAPK3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP013456HU  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P27361  |
| <b>Storage Buffer</b>    | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0  |
| <b>Product Type</b>      | Recombinant Proteins  |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | AAAAAQGGG GGEPRTTEGV GPGVPGEVEM VKGQPFVDVGP RYTQLQYIGE<br>GAYGMVSSAY DHVRKTRVAI KKISPFEHQY YCQRTLREIQ ILLRFRHENV<br>IGIRDILRAS TLEAMRDVYI VQDLMETDLY KLLKSQQLSN DHICYFLYQI<br>LRGLKYIHSA NVLHRDLKPS NLLINTTCDL KICDFGLARI ADPEHDHTGF<br>LTEYVATRWY RAPEIMLNSK GYTKSIDIWS VGCILAEMLS NRPIFPKGHY<br>LDQLNHILGI LGSPSQEDLN CIINMKARNY LQSLPSKTKV AWAKLFPKSD<br>SKALDLLDRM LTFNPNKRIT VEEALAHPYL EQYYDPTDEP VAEEPFTFAM<br>ELDDLPKERL KELIFQETAR FQPGVLEAP |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | MAPK3   |
| <b>Protein Names</b>     | Recommended name: Mitogen-activated protein kinase 3 Short name= MAP kinase 3 Short name= MAPK 3 EC= 2.7.11.24 Alternative name(s): ERT2 Extracellular signal-regulated kinase 1 Short name= ERK-1 Insulin-stimulat   |
| <b>Expression Region</b> | 2-332   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.