



Recombinant S- (hydroxymethyl)glutathione dehydrogenase (frmA)

Product Code	CSB-MP834778EOD
Abbreviation	frmA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X5J4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MKSRAAVAFA PGKPLEIVEI DVAPPKKGEV LIKVTHTGVC HTDAFTLSGD DPEGVFPVVL GHEGAGVVVE VGEGVTSVKP GDHVIPLYTA ECGECEFCRS GKTNLCVAVR ETQGKGLMPD GTTRFSYNGQ PLYHYMGCST FSEYTVVAEV SLAKINPEAN HEHVCLLGCG VTTGIGAVHN TAKVQPGDSV AVFGLGAIGL AVVQGARQAK AGRIIAIDTN PKKFDLARRF GATDCINPND YDKPIKDVL DINKWGIDHT FECIGNVNM RAALESAHRG WGQSVIIGVA GSGQEISTRP FQLVTGRVWK GSAFGGVKGR SQLPGMVEDA MKGDIDLEPF VTHTMLDEI NDAFDLMHEG KSIRTVIRY
Source	Mammalian cell
Target Names	frmA
Protein Names	Recommended name: S-(hydroxymethyl)glutathione dehydrogenase EC= 1.1.1.284 Alternative name(s): Alcohol dehydrogenase class-3 EC= 1.1.1.1 Alcohol dehydrogenase class-III Glutathione-dependent formaldehyde dehydrogenase Sh
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.