



Recombinant Uncharacterized protein ydaM (ydaM)

Product Code	CSB-EP834806EOD-B
Abbreviation	ydaM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X8Q3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	MITHNFNTLD LLTSPVWIVS PFEEQLIYAN SAARLLMQDL TFSQLRTGPY SVSSQKELPK YLSDLQHQHD IIEILTVQRK EEETALSCRL VLRELTETEP VIIFEGIEAP ATLGLKASRS ANYQRKKQGF YARFFLTNSA PMLLIDPSRD GQIVDANLAA LNFYGYNHET MCQKHTWEIN MLGRRVMPIM HEISHLPGGH KPLNFIHKLA DGSTRHVQTY AGPIEYGDK LMLCIVHDIT EQKRLEEQL HAAHHDAMTG LLNRRQFYHI TEPGQMQLA IAQDYSLLLI DTDRFKHIND LYGHSGKDEV LCALARTLES CARKGDLVFR WGGEEFVLLL PRTPLDTALS LAETIRVSA KVSISGLPRF TVSIGVAHHE GNESIDELFK RVDDALYRAK NDGRNRVLAA
Source	E.coli
Target Names	dgcM
Protein Names	Recommended name: Uncharacterized protein ydaM
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.