



Recombinant Vitamin B12-binding protein (btuF)

Product Code	CSB-BP834809EOD
Abbreviation	btuF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X8Z0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	APRVITLS PANTELAAFAA GITPVGVSYS SDYPPQAQKI EQVSTWQGMN LERIVALPKPD LVIWRGGNA ERQVDQLASL GIKVMWVDAT NIEQIANALR QLAPWSPQPD KAEQAAQSLL DQYAQLKAQY ADKPKKRVFL QFGINPPFTS GKESIQNQVL EVCGGENIFK DSRVPWPQVS REQVLARSPQ AIVITGGPDQ IPKIKQYWGE QLKIPVIPLT SDWFERASPR IILAAQQLCN ALSQVD
Source	Baculovirus
Target Names	btuF
Protein Names	Recommended name: Vitamin B12-binding protein
Expression Region	23-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.