



Recombinant L-aspartate oxidase (nadB)

Product Code	CSB-MP834818EOD
Abbreviation	nadB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8XA23
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MNTLLEHSCD VLIIGSGAAG LSLALRLADQ HQVIVLSKGP VTEGSTFYAQ GGIAAVFDET DSIDSHVEDT LIAGAGICDR HAVEFVASNA RSCVQWLIDQ GVLFDTHIQP NGEESYHLTR EGGHSHRRIL HAADATGREV ETTLVSKALN HPNIRVLERS NAVDLIISDK IGLPGTRRVV GAWVWNRNKE KVETCHAKAV VLATGGASKV YQYTTNPDIS SGDGIAMAWR AGCRVANLEF NQFHPTALYH PQARNFLLTE ALRGEAYLK RPDGTRFMPD FDVRGELAPR DIVARAI DHE MKRLGADCMF LDISHKPADF IRQHFPMIYE KLLGLGIDLT QEPVPIVPAA HYTCGGVMVD DHGRTDVEGL YAIGEVS YTG LHGANRMASN SLLECLVYGW SAAEDITRRM PYAHDISTLP PWDESRVENP DERVIIQHNW HELRLFMWDY VGIVRTTKRL ERALRRITML QQEIDEYYAH FRVSNLLEL RNLVQVAELI VRCAMMRKES RGLHFTLDYP ELLTHSGPSI LSPGNHYINR
Source	Mammalian cell
Target Names	nadB
Protein Names	Recommended name: L-aspartate oxidase Short name= LASPO EC= 1.4.3.16 Alternative name(s): Quinolate synthase B
Expression Region	1-540
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.