



Recombinant Cation efflux system protein CusB (cusB)

Product Code	CSB-EP834858EOD
Abbreviation	cusB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8XEH9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	EP PAEKTSTAER KILFWYDPMY PNTRFDKPGK SPFMDMDLVP KYADEESSAS GVRIDPTQTQ NLGVKTATVT RGPLTFAQSF PANVSYNEYQ YAIVQARAAG FIDKVYPLTV GDKVQKGAPL LDLTIPDWVE AQSEYLLRE TGGTATQTEG ILERLRLAGM PEADIRRLIA TQKIQTRFTL KAPIDGVITA FDLRAGMNIA KDNVVAKIQQ MDPVWVTAI PESIAWLVKD ASQFTLTVPA RPDKTLTIRK WTLLPGVDAA TRTLQLRLEV DNAMEALKPG MNAWLKLNTA SEPMLLIPSQ ALIDTGSEQR VITVDADGRF VPKRVAVFQA SQGV TALRSG LAEGEKVVSS GLFLIDSEAN ISGALERMRS ESATHAH
Source	E.coli
Target Names	cusB
Protein Names	Recommended name: Cation efflux system protein CusB
Expression Region	29-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.