



Recombinant Cyprinus carpio Dual specificity mitogen-activated protein kinase kinase 2 (map2k2)

Product Code	CSB-BP835446EQE
Abbreviation	map2k2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90321
Product Type	Recombinant Protein
Immunogen Species	Cyprinus carpio (Common carp)
Purity	≥85% (SDS-PAGE)
Sequence	MAPKRRPVPL IIAPTGEGQS TNIDAASEAN LEALQRKLGE LDLDEQQRKR LEAFLTQKAQ VGELKDEDFD PICELGAGNG GVVHKVRHKP SRLVMARKLI HLEIKPAIRN QIIRELQVLH ECNSPYIVGF YGAFYSDGEI SICMEHMDGG SLDQVLKEAR RIPEEILGKV SIAVLRGLVY LREKHQIMHR DVKPSNILVN SRGEIKLCDF GVSGQLIDSM ANSFVGTRSY MSPERLQGTH YSVQSDVWSM GLSLVELAIG RFPIPPDAK ELEAIFGRP V LDKGGAEGHS MSPRQRPPGR PVSGHGMDSR PAMAIFELLD YIVNEPPP K L PHGVFTTDFE EFVMKCLMKN PADRADLKML MGHTFIKRAE VEEVDFAGWM CKTMGLPQPS TPTTHSAE
Source	Baculovirus
Target Names	map2k2
Protein Names	Recommended name: Dual specificity mitogen-activated protein kinase kinase 2 Short name= MAP kinase kinase 2 Short name= MAPKK 2 EC= 2.7.12.2 Alternative name(s): ERK activator kinase 2 MAPK/ERK kinase 2 Short name=
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.