



Recombinant *Xenopus laevis* Suppressor of hairless protein homolog (rbpj), partial

Product Code	CSB-MP835525XBE
Abbreviation	rbpj
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91880
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MQPGIPKYTP SAIQLAPVVT GKFGERPQPK RLSTREAMRNY LKERGDQTVL ILHAKVAQKS YGNEKRFFCP PPCVYLMGSG WKKKKEQMER DGCSEQESQP CAFIGIGNSE QEMQQLNLEG KNYCTAKTLY ISDSDKRKH MLSVKMFYGN SDDIGVFLSK RIKVISKPSK KKQSLKNADL CIASGTKVAL FNRLRSQTVS TRYLHVEGGN FHASSQQWGA FYIHLLDDEE SEGEFTVRD GYIHYGQTVK LVCSVTGMAL PRLIIRKVDK QTALLDADDP VSQLHKCAFY LKDTERMYLC LSQERIIQFQ ATPCPKEPNK EMINDGASWT IISTDKAEYT FYEGMGPINA PVTPVPVVES LQLNGGGDVA MLELTGQNFT PNLRVWFGDV EAETMYRCAE SMLCVVPDIS AFREGWRWVR QPVQVPVTLV RNDGVIYSTS LTFTYTPEPG PRPHCSAAGA ILRANSSLLA SNEPNTNSEG SYTNISTNSA NVTSSSTAAVV S
Source	Mammalian cell
Target Names	rbpj
Protein Names	Recommended name: Suppressor of hairless protein homolog Alternative name(s): X-Su(H)
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.