



# Recombinant Human Phytanoyl-CoA hydroxylase-interacting protein (PHYHIP)

<b>Product Code</b>	CSB-MP835676HU
<b>Abbreviation</b>	PHYHIP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92561
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MELLSTPHSI EINNITCDSF RISWAMEDSD LERVTHYFID LNKKENKNSN KFKHRDVPTK LVAKAVPLPM TVRGHWFLSP RTEYSVAVQT AVKQSDGEYL VSGWSETVEF CTGDYAKEHL AQLQEKAEQI AGRMLRFSVF YRNHHKEYFQ HARTHCGNML QPYLKDNSGS HGSPTSGMLH GVFFSCNTEF NTGQPPQDSP YGRWRFQIPA QRLFNPSTNL YFADFYCMYT AYHYAILVLA PKGSLGDRFC RDRLPLLDIA CNKFLTCSVE DGELVFRHAQ DLILEIITYE PVDLSLGTLG EISGHQLMSL STADAKKDPS CKTCNISVGR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PHYHIP
<b>Protein Names</b>	Recommended name: Phytanoyl-CoA hydroxylase-interacting protein Alternative name(s): Phytanoyl-CoA hydroxylase-associated protein 1 Short name= PAHX-AP1 Short name= PAHXAP1
<b>Expression Region</b>	1-330
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.