



Recombinant Human Zinc finger protein neuro-d4 (DPF1)

Product Code	CSB-MP835697HU
Abbreviation	DPF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92782
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGGLSARPTA GRTDPAGTCW GQDPGSKMAT VIPGPLSLGE DFYREAIEHC RSYNARLCAE RSLRLPFLDS QTGVAQNNCY IWMEKTHRGP GLAPGQIYTY PARCWRKKRR LNILEDPRLR PCEYKIDCEA PLKKEGGLPE GPVLEALLCA ETGEKKIELK EEETIMDCQK QQLLEFP HDL EVEDLEDDIP RRKNRAKGKA YGIGGLRKRQ DTASLEDRDK PYVCDKFYKE LAWVPEAQRK HTAKKAPDGT VIPNGYCDFC LGGSKKTGCP EDLISCADCG RSGHP SCLQF TVNMTAAVRT YRWQCIECKS CSLCGTSEND GASWAGLTPQ DQLLFCDDCD RGYHMYCLSP PMAEPPEGSW SCHLCLRHLK EKASAYITLT
Source	Mammalian cell
Target Names	DPF1
Protein Names	Recommended name: Zinc finger protein neuro-d4 Alternative name(s): BRG1-associated factor 45B Short name= BAF45B D4, zinc and double PHD fingers family 1
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.