



Recombinant Arabidopsis thaliana 40S ribosomal protein S2-2 (RPS2B)

Product Code	CSB-MP835946DOA
Abbreviation	RPS2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93VB8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAERGGESGA ERGGDRGDFG RGFGGGRGGG RGRDRGPRGR GRRGGRASEE TKWVPVTKLG RLVADNKITK LEQIYLHSLP VKEYQIIDHL VGPTLKDEVM KIMPVQKQTR AGQRTRFKAF VVVGDNNGHV GLGVKCSKEV ATAIRGAIL AKLSVVPVRR GYWGNKIGKP HTVPCKVTGK CGSVTVRMVP APRGSIVAA RVPKKVLQFA GIDDVFTSSR GSTKTLGNFV KATFDCLQKT YGFLTPEFWK ETRFSRSPYQ EHTDFLSTKA VSATKVITEG EDQA
Source	Mammalian cell
Target Names	RPS2B
Protein Names	Recommended name: 40S ribosomal protein S2-2
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.