



Recombinant Arabidopsis thaliana tRNA (guanine (37)-N1)-methyltransferase 1 (At3g56120)

Product Code	CSB-EP835961DOA
Abbreviation	At3g56120
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93YU6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MFDESKFDVN LKLWALRIPR ELCKSASRIL NGYMLNMPRI KPITEDPTCE KTRLVILSES VKNADLSEIP EEKLNQLKKL SELEVPHSV TLGYSYWSAD HLLKQILPDG LDIPSSFETI GHIAHLNLHD ELLPFKDZIA KVIYDKNYPR IKTIVNKVGT ISNEFRVPKF EVLAGENGME TEVKQYGARF KLDYGLVYWN SRLEHEHMRL SSLFKPGETV CDMFAGIGPF AIPAAQKGCF VYANDLNPDS VRYLKINAKF NKVDDLICVH NMDARKFFSH LMAVSTCEDN LQSVADNDKT KEAAVSRGGE TNSSGEEIRE SNASINEPLG ANKKPSGTTK TENGVGKDCK SIEGHANKRL RQTLPIAKP WEHIDHVIMN LPASALQLD SFSNVIQKKY WKGPLPLIHC YCFIRASETT EFIIAEAETA LKFHIEDPVF HKVRDVAPNK AMFCLSFRLP EACLKQEE
Source	E.coli
Target Names	At3g56120
Protein Names	Recommended name: tRNA (guanine(37)-N1)-methyltransferase 1 EC=2.1.1.228 Alternative name(s): M1G-methyltransferase 1 tRNA [GM37] methyltransferase 1 tRNA methyltransferase 5 homolog 1
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.