



Recombinant Probable NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial (nduf-7), partial

Product Code	CSB-MP835979CXY
Abbreviation	nduf-7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94360
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	nduf-7
Protein Names	Recommended name: Probable NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-20kD Short name= CI-20kD NADH-ubiquinone oxidoreductase 20 kDa subunit
Expression Region	-
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.