



Recombinant Arabidopsis thaliana Peptidyl-prolyl cis-trans isomerase FKBP16-1, chloroplastic (FKBP16-1), partial

Product Code	CSB-MP835982DOA
Abbreviation	FKBP16-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q944B0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Source	Mammalian cell
Target Names	FKBP16-1
Protein Names	Recommended name: Peptidyl-prolyl cis-trans isomerase FKBP16-1, chloroplastic Short name= PPlase FKBP16-1 EC= 5.2.1.8 Alternative name(s): FK506-binding protein 16-1 Short name= AtFKBP16-1 Immunophilin FKBP16-1 Rotamas
Expression Region	-
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.