



Recombinant *Drosophila melanogaster* RNA-binding protein lark (lark)

| | |
|--------------------------|---|
| Product Code | CSB-BP835997DLU |
| Abbreviation | lark |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q94901 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila melanogaster</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MPGAGTFKLF IGNLDEKTQA TELRALFEKY GTVVECDVVK NYGFVHMETE QQGRDAIQNL NGYTLNEFAI KVEAAKSRRR PNTPTTKIFV GNLTDKTRAP EVRELFQKYG TVVECDIVRN YGFVHLCVVG DVQDAIKELN GRVVDGQPLK VQVSTSRVRP KPGMGDPEQC YRCGRSGHWS KECPRLYGSA GGGREPPSPL SAGGYRDRMY GRDPYPPPPP PPPFLRDRIM DGFRDYDYDD RRFEDSRDLY ERRYQTSRMR DFPPPPISRR EPMPPLPPTLS GSLRSCSVSR GYDTMFSRRS PPPPRSSNGM SRYGSPTPHG YEDFSRDAFD ERMISSRGMR GPSPPGRRYA PY |
| Source | Baculovirus |
| Target Names | lark |
| Protein Names | Recommended name: RNA-binding protein lark |
| Expression Region | 1-352 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |