



Recombinant *Drosophila melanogaster* RNA-binding protein lark (lark)

Product Code	CSB-EP835997DLU-B
Abbreviation	lark
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94901
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPGAGTFKLF IGNLDEKTQA TELRALFEKY GTVVECDVVK NYGFVHMETE QQGRDAIQNL NGYTLNEFAI KVEAAKSRRRA PNTPTTKIFV GNLTDKTRAP EVRELFQKYG TVVECDIVRN YGFVHLCVVG DVQDAIKELN GRVVDGQPLK VQVSTSRVVRP KPGMGDPEQC YRCGRSGHWS KECPRLYGSA GGGREPPSPL SAGGYRDRMY GRDPYPPPPP PPPFLRDRIM DGFRDYDYDD RRFEDSRDLY ERRYQTSRMR DFPPPPISRR EPMPPLPPTLS GSLRSCSVSR GYDTMFSRRS PPPPRSSNGM SRYGSPTPHG YEDFSRDAFD ERMISSRGMR GPSPPGRRYA PY
Source	E.coli
Target Names	lark
Protein Names	Recommended name: RNA-binding protein lark
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.