



Recombinant Glycine max Ferritin-2, chloroplastic

Product Code	CSB-MP836022GGV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94IC4
Product Type	Recombinant Protein
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Purity	>85% (SDS-PAGE)
Sequence	APAPLAGVI FEPFQELKGD YLAVPIAHNV XLARQNYADD SESAINEQIN VEYNVSYVYH ALFAYFDRDN IALKGLAKFF KESSEEREH AEQLIKYQNI RGGRVVLHPI TSPPEFEHS EKGDALYAME LALSLEKLTN EKLLHVHVA ERNNDPQXAD FIESEFLYEQ VKSIKIAEY VAQLRLVGKG HGVWHFDQKL LHDEDHV
Source	Mammalian cell
Protein Names	Recommended name: Ferritin-2, chloroplastic EC= 1.16.3.1 Alternative name(s): SFerH-2
Expression Region	52-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.