



Recombinant Putative T-box protein 7 (tbx-7)

Product Code	CSB-BP836116CXY
Abbreviation	tbx-7
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q95PX3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSEGYGRSLF SEKHFPHDRQ LLPPFTTTMS LGSFFSFPAS SAFFDRKKDD GVIDNPVEL VNRNLWSTFL ECGTEMIITK KGRRMFPLVK LKLSGLDKNS NYTIIMEMIS VDKLRYKFWN GNWIVAGVGE HHPLPTCFVH PQSPRSGEWW MTDGVDFKMA KLSNNPFNND GHIVLNSMHR YNPRFHIVRA DSSGQPILAS LKTFSFKETE FIAVTAYQND VVTCKIDNN PFAKGFRNID TVRKRKMNQI ISTPPASEDE EPEVKRTKSE IEAVMFSDPK VPQMSLQLIS QWQEALLSTI VQIPITNQKS SSERKSGFGV VDLLGSS
Source	Baculovirus
Target Names	tbx-7
Protein Names	Recommended name: Putative T-box protein 7
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>