



Recombinant Arabidopsis thaliana Peroxidase 9 (PER9)

Product Code	CSB-EP836151DOA
Abbreviation	PER9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96512
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	HPGLGFG WGSNSPIGGS FYSNLYPQFY QFSCPQADEI VMTVLEKAIA KEPRMAASLL RLHFHDCFVQ GCDASILLDD SATIRSEKNA GPNKNSVRGF QVIDEIKAKL EQACPQTVSC ADILALAARG STILSGGPSW ELPLGRRDSR TASLNGANTN IPAPNSTIQN LLTMFQRKGL NEEDLVSLSG GHTIGVARCT TFKQRLYNQN GNNQPDETLE RSYYYGLRSI CPPTGGDNNI SPLDLASPAR FDNTYFKLLL WGKGLLTSDE VLLTGNVGKT GALVKAYAED ERLFFQQFAK SMVNMGNIQP LTGFNGEIRK SCHVIN
Source	E.coli
Target Names	PER9
Protein Names	Recommended name: Peroxidase 9 Short name= Atperox P9 EC= 1.11.1.7 Alternative name(s): ATP18a
Expression Region	24-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.