



Recombinant Human Riboflavin kinase (RFK)

Product Code	CSB-EP836167HU
Abbreviation	RFK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q969G6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRHLPYFCRG QVVRGFGGRGS KQLGIPTANF PEQVVDNLPA DISTGIYYGW ASVGSQDVHK MVVSIGWNPY YKNTKKSMT HIMHTFKEDF YGEILNVAIV GYLRPEKNFD SLESLISAIQ GDIEEAKKRL ELPEHLKIKE DNFFQVSKSK IMNGH
Source	E.coli
Target Names	RFK
Protein Names	Recommended name: Riboflavin kinase EC= 2.7.1.26 Alternative name(s): ATP:riboflavin 5'-phosphotransferase Flavokinase
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.