



# Recombinant Human Zinc finger protein with KRAB and SCAN domains 4 (ZKSCAN4)

<b>Product Code</b>	CSB-BP836168HU
<b>Abbreviation</b>	ZKSCAN4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q969J2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAREPRKNAA LDAQSAEDQT GLLTVKVEKE EASALTAEVV APCSPARGPE RSRQRFRGFR YPEAAGPREA LSRLRELCGQ WLQPEMHSKE QILELLVLEQ FLTILPGNLQ SWVREQHPES GEEVVVLELEY LERQLDEPAP QVPVGDQGQE LLCCKMALLT QTQGSQSSQC QPMKALFKHE SLGSQPLHDR VLQVPGLAQQ GCCREDAMVA SRLTPGSQGL LKMEDVALTL TPGWTQLDSS QVNLYRDEKQ ENHSSLVSLG GEIQTCSRDL PPVKKLPEKE HGKICHLRED IAQIPTHAEA GEQEGRLQRK QKNAIGSRRH YCHECGKSFA QSSGLTKHRR IHTGEKPYEC EDCGKTFIGS SALVIHQRVH TGEKPYECEEE CGKVFSSHSSN LIKHQRTHTG EKPYECDDCG KTFQSQCSLL EHHKIHTGEK PYQCNMCGKA FRRNSHLLRH QRIHGDKNVQ NPEHGESWES QGRTESQWEN TEAPVSYKCN ECERSFTRNR SLIEHQKIHT GEKPYQCDCG GKGFTRTSYL VQHQRSHVVK KTLSSQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZKSCAN4
<b>Protein Names</b>	Recommended name: Zinc finger protein with KRAB and SCAN domains 4 Alternative name(s): P373c6.1 Zinc finger protein 307 Zinc finger protein 427
<b>Expression Region</b>	1-545
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.