



Recombinant Human Pentraxin-4 (PTX4)

Product Code	CSB-BP836179HU
Abbreviation	PTX4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96A99
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QEAAP VGPRKPFER LRRLEEQFRR FQEVTWTHLQ NIASNYNVSY NVDVRFRLA EESQAVAQAV NRSQASVQGE LAQLKAWVRK LQRRGRKVDT RLRALDTLG ERSQQRARER KAHKAQRDAL QDSLARLEGL VHSQGARLAA LEGRLPVAHP GTAALGPALV PTPTQPEELG PTSLKLQRDR QELRAASEHR GPPQDSSAPL QGRREPPASG SHRVLSTGAP KDPRQQAQWSP QVPGEICGVG PTLVFPNAST RNVVFLSPGF VTALRALSFC SWVRTASGRL GTLLSYATED NDNKLVHGR DSSLPGSIHF VIGDPAFREL PLQLLLDGQW HHICVIWTST QGRYWLHVDR RLVATGSRFR EGYEIPPGGS LVLGQEQDSV GGGFDSSEAF VGSMGLAIW DRALVPGEVA NLAIGKEFPT GAILTLANAA LAGGFVQGAN CTCLERCP
Source	Baculovirus
Target Names	PTX4
Protein Names	Recommended name: Pentraxin-4
Expression Region	26-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.