



Recombinant Human Clusterin-associated protein 1 (CLUAP1)

Product Code	CSB-YP836181HU
Abbreviation	CLUAP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96AJ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSFRDLRNFT EMMRALGYPR HISMENFRTP NFGVLVSEVLL WLVKRYEPQT DIPPDVDTEQ DRVFFIKAIA QFMATKAHIK LNTKKLYQAD GYAVKELLKI TSVLYNAMKT KGMEGSEIVE EDVNKFKFDL GSKIADLKAA RQLASEITSK GASLYDLLGM EVELREMRTE AIARPLEINE TEKVMRIAIAK EILTQVQKTK DLLNNVASDE ANLEAKIEKR KLELERNRKR LETLQSVRPC FMDEYEKTEE ELQKQYDTYL EKFNQLTYLE QQLEDHHRME QERFEEAKNT LCLIQNKLKE EEKRLLKSGS NDDSDIDIQE DDES DSELEE RRLPKPQTAM EMLMQGRP GK RIVGTMQGGD SDDNEDSEES EIDMEDDDDE DDDLEDESIS LSPTKPNRRV RKSEPLDESD NDF
Source	Yeast
Target Names	CLUAP1
Protein Names	Recommended name: Clusterin-associated protein 1
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.