



Recombinant Human Tripartite motif-containing protein 43 (TRIM43)

Product Code	CSB-MP836189HU
Abbreviation	TRIM43
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BQ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSDFSHAFQ KELTCVICLN YLVDPVTICC GHSFCRPCLC LSWEEAQSPA NCPACREPSP KMDFKTNILL KNLVTIARKA SLWQFLSSEK QICGTHRQTK KMFCDMDKSL LCLLCSNSQE HGAHKHHPPIE EAAEEHREKL LKQMRILWKK IQENQRNLYE EGRTAFLWRG NVVLRAQMIR NEYRKLHPVL HKEEKQHLE LNKEYQEIFQ QLQRSWVKMD QKSKHLKEMY QELMEMCHKP DVELLQDLGD IVARSESVLL HMPQPVNPEL TAGPITGLVY RLNRFRVEIS FHFEVTNHNI RLFEDVRSWM FRRGPLNSDR SDYFAAWGAR VFSFGKHYWE LDVDNSCDWA LGVCNNSWIR KNSTMVNSEID IFLLLCLKVD NHFNLLTTSP VFPHYIEKPL GRVGVFLDFE SGSVSFLNVT KSSLIWSYPA GSLTFPVRPF FYTGHR
Source	Mammalian cell
Target Names	TRIM43
Protein Names	Recommended name: Tripartite motif-containing protein 43
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.