



Recombinant Human 60S ribosomal export protein NMD3 (NMD3)

Product Code	CSB-MP836194HU
Abbreviation	NMD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96D46
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEYMAESTDR SPGHILCCEC GVPISPNPAN ICVACLRSKV DISQGIPKQV SISFCKQCQR YFQPPGTWIQ CALESRELLA LCLKKIKAPL SKVRLVDAGF VWTEPHSKRL KVKLTIQKEV MNGAILQQVF VVDYVVQSQM CGDCHRVEAK DFWKAVIQVR QKTLHKKTFY YLEQLILKYG MHQNTLRIKE IHDGLDFYYS SKQHAQKMVE FLQCTVPCRY KASQRLISQD IHSNTYNYKS TFSVEIVPIC KDNVVCLSPK LAQSLGNMNQ ICVCIRV TSA IHLIDPNTLQ VADIDGSTFW SHPFNSLCHP KQLEEFIVME CSIVQDIKRA AGAGMISKKH TLGEVWVQKT SEMNTDKQYF CRTHLGHLLN PGDLVLGFDL ANCNLNDEHV NKMNSDRVPD VVLIKKSYDR TKRQRRRNWK LKELAREREN MDTDDERQYQ DFLEDLEEDE AIRKNVNIYR DSAIPVESDT DDEGAPRISL AEMLEDLHIS QDATGEEGAS MLT
Source	Mammalian cell
Target Names	NMD3
Protein Names	Recommended name: 60S ribosomal export protein NMD3 Short name= hNMD3
Expression Region	1-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.