



Recombinant Human Zinc finger MYND domain-containing protein 19 (ZMYND19)

Product Code	CSB-MP836202HU
Abbreviation	ZMYND19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96E35
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTDFKLGIVR LGRVAGKTKY TLIDEQDIPL VESYSFEARM EVDADGNGAK IFAYAFDKNR GRGSGRLLHE LLWERHRGGV APGFQVVHLN AVTVDNRLDN LQLVPWGW RP KAEETSSKQR EQSLYWLAIQ QLPTDPIEEQ FPVLNVTRY NANGDVVEEE ENSCTYYECH YPPCTVIEKQ LREFNICGRC QVARYCGSQ CQQKDWPAHKK HCRERKRPFQ HELEPER
Source	Mammalian cell
Target Names	ZMYND19
Protein Names	Recommended name: Zinc finger MYND domain-containing protein 19 Alternative name(s): Melanin-concentrating hormone receptor 1-interacting zinc finger protein Short name= MCH-R1-interacting zinc finger protein
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.